



Ribbond®

bondable reinforcement ribbon



Instructions



Ribbond® bondable reinforcement

Ribbond-THM, Ribbond Ultra, Original Ribbond

Contents:

Ribbond starter kits contain:

Ribbond® bondable reinforcement ribbon, tinfoil and scissors.

Ribbond Scissors are only intended for cutting Ribbond fiber.

The scissors are not intended for other purposes.

Ribbond refill kits contain:

Ribbond® bondable reinforcement ribbon and tinfoil.

Statement of use: Ribbond reinforces dental resins.

Handling requirements: Before resin is applied, do not touch Ribbond with hands or latex gloves. Handle Ribbond with tweezers, pliers or cotton gloves. After resin is applied, handle Ribbond like the resin (instruments or powder free gloves).

Bonding Ribbond to composite: Wet Ribbond with unfilled bonding resin.
Imbed Ribbond within uncured filled composite.

Bonding Ribbond to acrylic: Wet Ribbond with a loose mix of liquid and powder.
Imbed Ribbond within uncured acrylic.

Polishing Ribbond: Fiber reinforcements do not polish. Cover Ribbond with filled composite or acrylic.

Disposal: No special considerations.

For Professional Use Only

Ribbond® Bondable Reinforcement Ribbon | Made in the USA

US Patent - 9138298 Europe - 2731537

Germany - 6020120439866 - Japan - 特許第5922770号



Ribbond, Inc.
2200 South St., Racine WI 53404 USA
800-624-4554
ribbond@ribbond.com
www.ribbond.com



European Representative:
Sigma Dental Systems
Heideland 22
24976 Handewitt
Germany



Ribbond®

bondable reinforcement ribbon

CATALOG NO.	DESCRIPTION	CONTENTS
Ribbond Bondable Reinforcement Ribbon		
AST	Assorted Starter Kit	One scissors, one tin foil, and one 3 x 22 cm Ribbond.
CUST	Custom Assorted Starter Kit	One scissors, one tin foil and one 1 x 68 cm Ribbond.
KIT2	2 mm Starter Kit	One scissors, one tin foil, and one 3 x 22 cm Ribbond.
KIT3	3 mm Starter Kit	One scissors, one tin foil and one 1 x 68 cm Ribbond.
KIT4	4 mm Starter Kit	One scissors, one tin foil, and one 3 x 22 cm Ribbond.
KITO	Orthodontic Starter Kit	One scissors, one tin foil and one 1 x 68 cm Ribbond.
REAST	Assorted Refill	One tin foil and one 3 x 22 cm Ribbond.
CUSTREAST	Custom Assorted Refill	One tin foil and one 1 x 68 cm Ribbond.
RE2	2 mm Refill	One tin foil and one 3 x 22 cm Ribbond.
RE3	3 mm Refill	One tin foil and one 1 x 68 cm Ribbond.
RE4	4 mm Refill	One tin foil and one 3 x 22 cm Ribbond.
RE9	9 mm Refill	One tin foil and one 1 x 45 cm Ribbond.
Ribbond THM (Thinner Higher Modulus)		
ASTTHM	Assorted Starter Kit - THM	One scissors, one tin foil, and one 3 x 22 cm Ribbond THM.
CUSTTHM	Custom Starter Kit - THM	One scissors, one tin foil and one 1 x 68 cm Ribbond THM
KIT2THM	2 mm Starter Kit - THM	One scissors, one tin foil, and one 3 x 22 cm Ribbond THM.
KIT3THM	3 mm Starter Kit - THM	One scissors, one tin foil and one 1 x 68 cm Ribbond THM.
KIT4THM	4 mm Starter Kit - THM	One scissors, one tin foil, and one 3 x 22 cm Ribbond THM.
KITOTHM	Orthodontic Starter Kit -THM	One scissors, one tin foil and one 1 x 68 cm Ribbond THM.
REASTTHM	Assorted Refill -THM	One tin foil and one 3 x 22 cm Ribbond THM.
CUSTRETHM	Custom Assorted Refill -THM	One tin foil and one 1 x 68 cm Ribbond THM.
RE2THM	2 mm Refill -THM	One tin foil and one 3 x 22 cm Ribbond THM.
RE3THM	3 mm Refill -THM	One tin foil and one 1 x 68 cm Ribbond THM.
RE4THM	4 mm Refill -THM	One tin foil and one 3 x 22 cm Ribbond THM.
REOTHM	Ortho Refill -THM	One tin foil and one 1 x 68 cm Ribbond THM.
Ribbond Triaxial		
ASTTRI	Assorted Triaxial Starter Kit	One scissors, one tin foil and one 4 x 15 cm Ribbond Triaxial.
REASTTRI	Assorted Triaxial Refill	One tin foil and one 4 x 15 cm Ribbond Triaxial.